


Anti-Adra2A Rabbit Polyclonal Antibody For Wb - P08913

Item Number: CM0025251



Introduction

Anti-ADRA2A rabbit polyclonal antibody, validated for WB and ELISA. Detects a 51kDa band and cross-reacts with human, mouse, and rat ADRA2A. Ideal for GPCR signaling, cardiovascular, and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	ADRA2A Rabbit pAb
Remarks/Alias	ADRA2; ADRAR; ZNF32; ADRA2R; ALPHA2AAR; ADRA2A
Species	Human
GeneID (human)	150
GeneID	150
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 240-380 of human ADRA2A (NP_000672.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08913
Protein Weight	51kDa
Shipping	Ice bag