


Anti-Adra2B Polyclonal Antibody For Wb And Elisa - P18089

Item Number: CM0025643



Introduction

Rabbit polyclonal antibody against human Alpha-2B adrenergic receptor (ADRA2B). Validated for Western blot and ELISA; reacts with human, mouse, rat. Ideal for studying blood pressure regulation, GPCR signaling, and metabolic processes.

[Learn More](#)

Parameter	Value
Product Name	ADRA2B Rabbit pAb
Remarks/Alias	FAME2; ADRA2L1; ADRARL1; ADRA2RL1; ALPHA2BAR; alpha-2BAR; ADRA2B
Species	Human
GeneID (Human)	151
GeneID	151
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ADRA2B (NP_000673.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18089
Protein Weight	50kDa
Shipping	Ice bag