


Anti-Adra2C Rabbit Polyclonal Antibody For Wb, Elisa - P18825

Item Number: CM0021483



Introduction

Rabbit polyclonal antibody targeting human ADRA2C (Alpha-2C adrenergic receptor). Validated for WB, ELISA; cross-reacts with human and mouse. Suitable for cardiovascular and neurological research on blood pressure, lipolysis, insulin release, and neurotransmission.

[Learn More](#)

Parameter	Value
Product Name	ADRA2C Rabbit pAb
Remarks/Alias	ADRA2L2; ADRARL2; ADRA2RL2; ALPHA2CAR; ADRA2C
Species	Human
GeneID (human)	152
GeneID	152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-380 of human ADRA2C (NP_000674.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P18825
Protein Weight	50kDa
Shipping	Ice bag