


Anti-Adrb2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P07550

Item Number: CM0019647



Introduction

High-quality anti-ADRB2 rabbit polyclonal antibody validated in WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying beta-2 adrenergic receptor signaling in cardiovascular and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	ADRB2 Rabbit pAb
Remarks/Alias	BAR; B2AR; ADRBR; ADRB2R; BETA2AR; ADRB2
Species	Human
GeneID (human)	154
GeneID	154
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ADRB2 (NP_000015.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07550
Protein Weight	46kDa
Shipping	Ice bag