


Anti-Adrb3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13945

Item Number: CM0027361



Introduction

High-quality rabbit polyclonal antibody against human ADRB3 (Beta-3 adrenergic receptor). Validated for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Ideal for thermogenesis, lipolysis, and GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	ADRB3 Rabbit pAb
Remarks/Alias	BETA3AR; ADRB3
Species	Human
Gene ID (Human)	155
Gene ID	155
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ADRB3 (NP_000016.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13945
Protein Weight	44kDa
Shipping	Ice bag