

Anti-Adrbk2 (Grk3, Bark2) Polyclonal Antibody For Wb, Ihc-P, Elisa - P35626

Item Number: CM0023564

Introduction

High-quality rabbit polyclonal antibody against human ADRBK2 (GRK3/BARK2) for WB, IHC-P, ELISA. Validated in human, mouse, rat. Investigate GPCR signaling and desensitization.

[Learn More](#)

Parameter	Details
Product Name	ADRBK2 Rabbit pAb
Remarks/Alias	BARK2; ADRBK2
Species	Human
GeneID (human)	157
GeneID	157
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ADRBK2 (NP_005151.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35626
Protein Weight	80kDa
Shipping	Ice bag