


Anti-Adrm1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16186

Item Number: CM0030215



Introduction

Rabbit polyclonal antibody targeting ADRM1 (proteasomal ubiquitin receptor, Rpn13). Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat. Supports proteasome and ubiquitin research.

[Learn More](#)

Parameter	Value
Product Name	ADRM1 Rabbit pAb
Remarks/Alias	ARM1; ARM-1; GP110; PSMD16; ADRM1
Species	Human
GeneID (Human)	11047
GeneID	11047
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 108-407 of human ADRM1 (NP_008933.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16186
Protein Weight	42kDa
Shipping	Ice bag