


Anti-Adsl Polyclonal Antibody For Wb, If/Icc, Elisa - P30566

Item Number: CM0025815



Introduction

Anti-ADSL polyclonal antibody validated for WB, IF/ICC, ELISA. Cross-reacts with Human, Mouse, Rat. Targets purine biosynthesis enzyme adenylosuccinate lyase. Ideal for metabolic research.

[Learn More](#)

Parameter	Value
Product Name	ADSL Rabbit pAb
Remarks/Alias	ASL; AMPS; ASASE; ADSL
Species	Human
GeneID (human)	158
GeneID	158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human ADSL (NP_000017.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30566
Protein Weight	55kDa
Shipping	Ice bag