

Anti-Adss 2 Rabbit Monoclonal Antibody For Wb And Elisa - P30520

Item Number: CM0002377



Introduction

Rabbit monoclonal antibody targeting human adenylosuccinate synthetase isozyme 2 (AdSS 2), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for purine metabolism research.

[Learn More](#)

Parameter	Value
Product Name	AdSS 2 Rabbit mAb
Remarks/Alias	ADEH; ADSS; ADSS 2; AdSS 2
Species	Human
GeneID (Human)	159
GeneID	159
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30520
Protein Weight	50kDa
Shipping	Ice bag