

Anti-Aebp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8lux7

Item Number: CM0019563

Introduction

High-quality AEBP1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for extracellular matrix research, adipocyte biology, and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	AEBP1 Rabbit pAb
Remarks/Alias	ACLP; AEBP1
Species	Human
GeneID (Human)	165
GeneID	165
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 999-1158 of human AEBP1 (NP_001120.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IUX7
Protein Weight	131kDa
Shipping	Ice bag