

Anti-Aebp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Zn18

Item Number: CM0030739



Introduction

Rabbit polyclonal antibody targeting human AEBP2 (Zinc finger protein AEBP2), validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. AEBP2 is an accessory subunit of PRC2 that mediates H3K27 trimethylation and transcriptional repression. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	AEBP2 Rabbit pAb
Remarks/Alias	AEBP2
Species	Human
GeneID (Human)	121536
GeneID	121536
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-503 of human AEBP2 (NP_694939.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZN18
Protein Weight	54kDa
Shipping	Ice bag