

Anti-Aer61/Eogt Rabbit Monoclonal Antibody For Wb And Elisa - Q5NDL2

Item Number: CM0003734

Introduction

High-quality rabbit monoclonal antibody against human EGF domain-specific O-linked N-acetylglucosamine transferase (EOGT/AER61). Validated for WB and ELISA, cross-reacts with human, mouse, rat. Ideal for glycobiology and protein modification research.

[Learn More](#)

Parameter	Value
Product Name	AER61 Rabbit mAb
Remarks/Alias	AOS4; AER61; EOGT1; C3orf64
Species	Human
GeneID (Human)	285203
GeneID	285203
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human AER61 (NP_775925.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5NDL2
Protein Weight	62kDa
Shipping	Ice bag