


Anti-Aff4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uhb7

Item Number: CM0028559



Introduction

Rabbit polyclonal anti-AFF4 antibody (SWISS: Q9UHB7) validated for WB and ELISA, suitable for transcription elongation complex studies in human and mouse samples. Ideal for SEC complex research.

[Learn More](#)

Parameter	Value
Product Name	AFF4 Rabbit pAb
Remarks/Alias	MCEF; CHOPS; AF5Q31; AFF4
Species	Human
Gene ID (Human)	27125
Gene ID	27125
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human AFF4 (NP_055238.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UHB7
Protein Weight	127kDa
Shipping	Ice bag