


Anti-Afg3L2 Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q9Y4W6

Item Number: CM0010031



Introduction

Rabbit monoclonal antibody against human AFG3L2, validated for WB, IP, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial protease and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	AFG3L2 Rabbit PolymAb®
Remarks/Alias	OPA12; SCA28; SPAX5
Species	Human
Gene ID (Human)	10939
Gene ID	10939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 538-797 of human AFG3L2 (NP_006787.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4W6
Protein Weight	89kDa
Shipping	Ice bag