


Anti-Afg3L2 (Opa12/Sca28/Spax5) Monoclonal Antibody For Wb, Ihc-P, If/Icc - Q9Y4W6

Item Number: CM0009080



Introduction

Anti-AFG3L2 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, and ELISA. Reacts with human, mouse, rat AFG3L2 (89 kDa), a mitochondrial m-AAA protease implicated in spinocerebellar ataxia.

[Learn More](#)

Parameter	Value
Product Name	AFG3L2 Rabbit mAb
Remarks/Alias	OPA12; SCA28; SPAX5
Species	Human
GeneID (Human)	10939
GeneID	10939
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 538-797 of human AFG3L2 NP_006787.2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4W6
Protein Weight	89kDa
Shipping	Ice bag