


Anti-Afg3L2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Y4W6

Item Number: CM0022339



Introduction

Rabbit anti-AFG3L2 polyclonal antibody validated for WB, IHC-P, IP, ELISA; cross-reacts with human, mouse, rat. Recombinant immunogen. Ideal for mitochondrial inner membrane m-AAA protease studies and spinocerebellar ataxia research.

[Learn More](#)

Parameter	Value
Product Name	AFG3L2 Rabbit pAb
Remarks/Alias	OPA12; SCA28; SPAX5; AFG3L2
Species	Human
Gene ID (Human)	10939
Gene ID	10939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 538-797 of human AFG3L2 (NP_006787.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4W6
Protein Weight	89kDa
Shipping	Ice bag