


Anti-Afg3L2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Y4W6

Item Number: CM0021149



Introduction

Rabbit polyclonal antibody against human AFG3L2 (m-AAA protease), validated for WB, IHC-P, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial quality control and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	AFG3L2 Rabbit pAb
Remarks/Alias	OPA12; SCA28; SPAX5; AFG3L2
Species	Human
GeneID (Human)	10939
GeneID	10939
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 538-797 of human AFG3L2 (NP_006787.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4W6
Protein Weight	89kDa
Shipping	Ice bag