

Anti-Afm Rabbit Polyclonal Antibody For Wb, Elisa - P43652

Item Number: CM0017198



Introduction

Rabbit polyclonal antibody against human afamin (AFM), validated for WB and ELISA with cross-reactivity to rat. Ideal for Wnt signaling and vitamin E transport studies.

[Learn More](#)

Parameter	Value
Product Name	AFM Rabbit pAb
Remarks/Alias	ALF; ALB2; ALBA; AFM
Species	Human
GeneID (Human)	173
GeneID	173
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-599 of human AFM (NP_001124.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P43652
Protein Weight	69kDa
Shipping	Ice bag