

Anti-Ager Monoclonal Antibody For Wb - Q15109

Item Number: CM0008675



Introduction

Rabbit recombinant anti-AGER (RAGE) monoclonal antibody validated for WB, IF, IHC, and ELISA. Reacts with human, mouse, rat. Ideal for RAGE signaling, inflammation, and age-related disease studies.

[Learn More](#)

Parameter	Value
Product Name	AGER Rabbit PolymAb®
Remarks/Alias	RAGE; sRAGE; SCARJ1
Species	Human
Gene ID (Human)	177
Gene ID	177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-344 of human AGER. (NP_001127.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15109
Protein Weight	43kDa
Shipping	Ice bag