

Anti-Ager Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q15109

Item Number: CM0006575



Introduction

Rabbit monoclonal antibody specific for human AGER/RAGE (43kDa). Suitable for WB, IF-P, and ELISA. Ideal for studying inflammation, diabetes complications, and receptor-mediated signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	AGER Rabbit mAb
Remarks/Alias	RAGE; sRAGE; SCARJ1; AGER
Species	Human
GeneID	177
Immunogen	Recombinant protein of human AGER
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	Q15109
Protein Weight	43kDa
Shipping	Ice bag