

Anti-Ager (Rage) Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q15109

Item Number: CM0006432



Introduction

Recombinant rabbit monoclonal antibody against human AGER (RAGE). Validated for WB, IF-P, IHC-P, ELISA applications. Cross-reacts with mouse and rat. Ideal for inflammation, diabetes, and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	AGER Rabbit mAb
Remarks/Alias	RAGE; sRAGE; SCARJ1; AGER
Species	Human
GeneID	177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-344 of human AGER (NP_001127.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15109
Protein Weight	43kDa
Shipping	Ice bag