

Anti-Ager (Rage) Polyclonal Antibody For Wb, Elisa - Q15109

Item Number: CM0016419



Introduction

Rabbit polyclonal antibody against human AGER (RAGE), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Targets the extracellular domain of RAGE, a pattern recognition receptor implicated in inflammation, diabetes complications, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	AGER Rabbit pAb
Remarks/Alias	RAGE; sRAGE; SCARJ1; AGER
Species	Human
GeneID	177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-210 of human AGER (NP_001127.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15109
Protein Weight	43kDa
Shipping	Ice bag