

Anti-Ager (Rage) Rabbit Polyclonal Antibody For Wb And Elisa - Q15109

Item Number: CM0016131



Introduction

Rabbit polyclonal antibody against AGER (RAGE) for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for research on advanced glycation end products, diabetic complications, and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	AGER Rabbit pAb
Remarks/Alias	RAGE; sRAGE; SCARJ1; AGER
Species	Human
GeneID (Human)	177
GeneID	177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-330 of human AGER (NP_001127.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15109
Protein Weight	43kDa
Shipping	Ice bag