


Anti-Agfg1 Polyclonal Antibody For Wb And Elisa - P52594

Item Number: CM0030791



Introduction

Rabbit polyclonal antibody targeting human AGFG1, validated for Western blot and ELISA. Suitable for investigating HIV-1 Rev-dependent RNA export, acrosome biogenesis, and RNA trafficking. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	AGFG1 Rabbit pAb
Remarks/Alias	HRB, RAB, RIP, AGFG1
Species	Human
GeneID (Human)	3267
GeneID	3267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human AGFG1 (NP_004495.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52594
Protein Weight	58kDa
Shipping	Ice bag