

Anti-Aggrecan Rabbit Polyclonal Antibody For Wb, Ihc, If, Elisa - P16112

Item Number: CM0027279



Introduction

Rabbit polyclonal antibody against human Aggrecan, validated for WB, IF-P, IHC-P, and ELISA. Detects Aggrecan in human, mouse, and rat samples. Ideal for cartilage matrix research.

[Learn More](#)

Parameter	Value
Product Name	Aggrecan Rabbit pAb
Remarks/Alias	AGC1; SEDK; AGCAN; CSPG1; MSK16; CSPGCP; SSOAOD; Aggrecan
Species	Human
GeneID (Human)	176
GeneID	176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 420-575 of human Aggrecan (NP_037359.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16112
Protein Weight	261kDa
Shipping	Ice bag