

Anti-Agk Mouse Monoclonal Antibody For Wb And Elisa - Q53H12

Item Number: CM0002497



Introduction

High-quality anti-AGK mouse monoclonal antibody validated for WB and ELISA applications. Specific for human acylglycerol kinase, a mitochondrial lipid kinase involved in LPA/PA production and protein import.

[Learn More](#)

Parameter	Value
Product Name	AGK Mouse mAb
Remarks/Alias	MULK; CATC5; CTRCT38; MTDPS10; AGK
Species	Human
GeneID (Human)	55750
GeneID	55750
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-422 of human AGK (NP_060708.1)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q53H12
Protein Weight	47kDa
Shipping	Ice bag