

Anti-Agk Rabbit Polyclonal Antibody For Wb, Elisa - Q53H12

Item Number: CM0030337



Introduction

Rabbit polyclonal antibody against human AGK, suitable for WB and ELISA. Reacts with human and rat. AGK is a mitochondrial lipid kinase involved in LPA/PA production and TIM22 complex-mediated protein import. Ideal for metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	AGK Rabbit pAb
Remarks/Alias	MULK; CATC5; CTRCT38; MTDPS10; AGK
Species	Human
Gene ID (Human)	55750
Gene ID	55750
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q53H12
Protein Weight	47kDa
Shipping	Ice bag