


Anti-Agl (Gde) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P35573

Item Number: CM0004660



Introduction

This Anti-AGL (GDE) Rabbit Monoclonal Antibody is validated for Western blot, immunofluorescence, and ELISA. It cross-reacts with human, mouse, and rat samples, supporting glycogen debranching enzyme function and localization studies.

[Learn More](#)

Parameter	Value
Product Name	AGL Rabbit mAb
Remarks/Alias	GDE; AGL
Species	Human
GeneID (Human)	178
GeneID	178
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human AGL (P35573).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35573
Protein Weight	175kDa
Shipping	Ice bag