

Anti-Agmat Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bse5

Item Number: CM0007102

Introduction

Rabbit monoclonal antibody against human AGMAT, suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Detects mitochondrial guanidino acid hydrolase, ~38 kDa.

[Learn More](#)

Parameter	Value
Product Name	AGMAT Rabbit mAb
Remarks/Alias	AGMAT
Species	Human
Gene ID	79814
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-352 of human AGMAT (NP_079034.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BSE5
Protein Weight	38 kDa
Shipping	Ice bag