

Anti-Agmat Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bse5

Item Number: CM0013145

Introduction

Anti-AGMAT rabbit polyclonal antibody validated for WB, IF/ICC, ELISA applications with human and rat species reactivity. Recognizes endogenous AGMAT (guanidino acid hydrolase), a 38 kDa mitochondrial enzyme in arginine metabolism and cellular detoxification.

[Learn More](#)

Parameter	Value
Product Name	AGMAT Rabbit pAb
Remarks/Alias	AGMAT
Species	Human
Gene ID	79814
Immunogen	Recombinant fusion protein corresponding to amino acids 1-352 of human AGMAT (NP_079034.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BSE5
Protein Weight	38kDa
Shipping	Ice bag