

Anti-Ago1 Polyclonal Antibody For Wb, Elisa - Q9UL18

Item Number: CM0024837



Introduction

AGO1 Rabbit pAb is a rabbit polyclonal antibody raised against human AGO1 N-terminus. Validated in WB and ELISA, it detects AGO1 in human, mouse, and rat samples, making it suitable for RNAi pathway studies.

[Learn More](#)

Parameter	Value
Product Name	AGO1 Rabbit pAb
Remarks/Alias	Q99; EIF2C; hAgo1; EIF2C1; GERP95; NEDLBAS; AGO1
Species	Human
GeneID (Human)	26523
GeneID	26523
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human AGO1 (NP_036331.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UL18
Protein Weight	97kDa
Shipping	Ice bag