

Anti-Agpat1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q99943

Item Number: CM0028859

Introduction

AGPAT1 Rabbit Polyclonal Antibody validated for WB, IP, and ELISA. Detects human, mouse, and rat AGPAT1, a key acyltransferase in phospholipid biosynthesis. Suitable for lipid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	AGPAT1 Rabbit pAb
Remarks/Alias	G15; LPAATA; LPLAT1; 1-AGPAT1; LPAAT-alpha; AGPAT1
Species	Human
GeneID	10554
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 204-283 of human AGPAT1 (NP_006402.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99943
Protein Weight	32kDa
Shipping	Ice bag