

Anti-Agpat2 Polyclonal Antibody For Wb, If/Icc, Elisa - O15120

Item Number: CM0030364

Introduction

Rabbit polyclonal antibody targeting human AGPAT2 (1-acylglycerol-3-phosphate O-acyltransferase beta), validated for WB, IF/ICC, and ELISA. Reacts with human and mouse. Suitable for lipid metabolism and lipodystrophy research.

[Learn More](#)

Parameter	Specification
Product Name	AGPAT2 Rabbit pAb
Remarks/Alias	BSCL; BSCL1; LPAAB; LPLAT2; 1-AGPAT2; LPAAT-beta; AGPAT2
Species	Human
GeneID (human)	10555
GeneID	10555
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 199-278 of human AGPAT2 (NP_006403.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15120
Protein Weight	31kDa
Shipping	Ice bag