

Anti-Agpat4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nrz5

Item Number: CM0021319

Introduction

Rabbit polyclonal antibody against human AGPAT4, validated for WB and ELISA. Detects mouse. Suitable for studies of phospholipid biosynthesis, lysophosphatidic acid acyltransferase activity, and lipid metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	AGPAT4 Rabbit pAb
Remarks/Alias	LPAAT4; LPLAT4; 1-AGPAT4; dJ473J16.2; LPAAT-delta; AGPAT4
Species	Human
GeneID (human)	56895
GeneID	56895
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 270-378 of human AGPAT4 (NP_064518.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NRZ5
Protein Weight	44kDa
Shipping	Ice bag