

Anti-Agps Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 000116

Item Number: CM0017742

Introduction

Rabbit polyclonal antibody targeting human AGPS, suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse. For peroxisomal ether lipid biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	AGPS Rabbit pAb
Remarks/Alias	ADAS, ADPS, RCDP3, ADAP-S, ADHAPS, ALDHPSY, AGPS
Species	Human
Gene ID (Human)	8540
Gene ID	8540
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 379-658 of human AGPS (NP_003650.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00116
Protein Weight	73kDa
Shipping	Ice bag