

Anti-Agr2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O95994

Item Number: CM0003754



Introduction

High-quality AGR2 rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Detects human and mouse AGR2, a proto-oncogene involved in mucus secretion and cell adhesion. Ideal for cancer and secretory pathway research.

[Learn More](#)

Parameter	Specification
Product Name	AGR2 Rabbit mAb
Remarks/Alias	AG2; AG-2; HPC8; GOB-4; HAG-2; RIFTD; XAG-2; PDIA17; HEL-S-116; AGR2
Species	Human
GeneID (human)	10551
GeneID	10551
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human AGR2 (O95994).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95994
Protein Weight	20kDa
Shipping	Ice bag