

Anti-Agtr2 Rabbit Monoclonal Antibody For Wb, Elisa - P50052

Item Number: CM0011544



Introduction

Anti-AGTR2 Rabbit Monoclonal Antibody (WB, ELISA) raised against a synthetic peptide from human AGTR2. Cross-reacts with human, mouse, and rat. For angiotensin II receptor research.

[Learn More](#)

Parameter	Value
Product Name	AGTR2 Rabbit mAb
Remarks/Alias	AT2; ATGR2; MRX88; AGTR2
Species	Human
Gene ID	186
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 264-363 of human AGTR2 (P50052).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50052
Protein Weight	41kDa
Shipping	Ice bag