

Anti-Agtr2 Rabbit Polyclonal Antibody For Wb And Elisa - P50052

Item Number: CM0021179



Introduction

High-quality rabbit polyclonal antibody targeting human AGTR2 (Angiotensin II type-2 receptor). Validated for Western blot and ELISA. Suitable for research on vasoconstriction, cell proliferation inhibition, and GPCR signaling.

[Learn More](#)

Parameter	Value
Product Name	AGTR2 Rabbit pAb
Remarks/Alias	AT2, ATGR2, MRX88, AGTR2
Species	Human
Gene ID (Human)	186
Gene ID	186
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human AGTR2 (NP_000677.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P50052
Protein Weight	41kDa
Shipping	Ice bag