

Anti-Agtrap Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Rw13

Item Number: CM0016489

Introduction

Anti-AGTRAP rabbit polyclonal antibody validated for WB, IHC-P, and ELISA in human, mouse, and rat. Targets the 17 kDa type-1 angiotensin II receptor-associated protein, a negative regulator of angiotensin II signaling. Ideal for cardiovascular and renal disease research.

[Learn More](#)

Parameter	Value
Product Name	AGTRAP Rabbit pAb
Remarks/Alias	ATRAP; AGTRAP
Species	Human
GeneID (Human)	57085
GeneID	57085
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 93-152 of human AGTRAP (NP_065083.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6RW13
Protein Weight	17 kDa
Shipping	Ice bag