


Anti-Aha1 Monoclonal Antibody For Wb, Elisa - 095433

Item Number: CM0008673



Introduction

High-quality rabbit monoclonal antibody targeting human AHA1/p38 (activator of HSP90 ATPase). Validated for Western blot and ELISA with cross-reactivity to mouse and rat. Ideal for chaperone and stress response studies.

[Learn More](#)

Parameter	Value
Product Name	AHA1 Rabbit mAb
Remarks/Alias	p38, AHA1, hAha1, C14orf3
Species	Human
Gene ID (Human)	10598
Gene ID	10598
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 239-338 of human AHA1 (O95433).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95433
Protein Weight	38kDa
Shipping	Ice bag