


Anti-Ahcy Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P23526

Item Number: CM0029359



Introduction

AHCY Rabbit Polyclonal Antibody (SAHH) validated for WB, IF/ICC, ELISA against human, mouse, rat. Targets adenosylhomocysteinase, 48 kDa, essential for one-carbon metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	AHCY Rabbit pAb
Remarks/Alias	SAHH; adoHcyase; AHCY
Species	Human
GeneID (human)	191
GeneID	191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 253-432 of human AHCY (NP_000678.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23526
Protein Weight	48kDa
Shipping	Ice bag