


Anti-Ahcy1 Rabbit Polyclonal Antibody For Wb, Elisa - O43865

Item Number: CM0030441



Introduction

Rabbit polyclonal antibody against AH CYL1 (IRBIT) for Western blot and ELISA. Detects human, mouse, and rat S-adenosylhomocysteine hydrolase-like protein 1, involved in calcium signaling, fluid secretion, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	AH CYL1 Rabbit pAb
Remarks/Alias	DCAL; IRBIT; PPP1R78; PRO0233; XPVKONA; AH CYL1
Species	Human
GeneID (Human)	10768
GeneID	10768
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-260 of human AH CYL1 (NP_006612.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43865
Protein Weight	59kDa
Shipping	Ice bag