

Anti-Ahr Rabbit Monoclonal Antibody For Wb, Elisa - P35869

Item Number: CM0009038



Introduction

High-quality AHR Rabbit Monoclonal Antibody for WB and ELISA. Targets human aryl hydrocarbon receptor (bHLHe76, P35869). Ideal for xenobiotic response, cancer, and circadian rhythm research.

[Learn More](#)

Parameter	Value
Product Name	AHR Rabbit mAb
Remarks/Alias	RP85; bHLHe76; AHR
Species	Human
GeneID (Human)	196
GeneID	196
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 619-848 of human AHR (NP_001612.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35869
Protein Weight	96kDa
Shipping	Ice bag