

Anti-Ahr Rabbit Polyclonal Antibody For Wb And Elisa - P30561

Item Number: CM0029345



Introduction

Anti-AHR Rabbit Polyclonal Antibody tested for Western blotting and ELISA. Recognizes mouse AHR, a ligand-activated transcription factor involved in xenobiotic metabolism, immunity, and circadian rhythms. Ideal for cancer and toxicology research.

[Learn More](#)

Parameter	Value
Product Name	AHR Rabbit pAb
Remarks/Alias	Ah; In; Ahh; Ahre; bHLHe76
Species	Mouse
GeneID (Human)	196
GeneID	11622
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 706-805 of mouse AHR (NP_038492.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P30561
Protein Weight	95 kDa
Shipping	Ice bag