

Anti-Ahr Rabbit Polymab® Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35869

Item Number: CM0007128



Introduction

Rabbit monoclonal antibody against human aryl hydrocarbon receptor (AHR, bHLHe76). Suitable for WB, IHC-P, IF/ICC, ELISA. Ideal for studies on xenobiotic response, circadian rhythm, and cancer. Recognizes a 96 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	AHR Rabbit PolymAb®
Remarks/Alias	RP85; bHLHe76
Species	Human
Gene ID (Human)	196
Gene ID	196
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 619-848 of human AHR Rabbit PolymAb® (NP_001612.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P35869
Protein Weight	96kDa
Shipping	Ice bag