


Anti-Ahsp Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nzd4

Item Number: CM0026855



Introduction

Rabbit polyclonal antibody against human AHSP, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for erythropoiesis, hemoglobin stability, and beta-thalassemia research.

[Learn More](#)

Parameter	Value
Product Name	AHSP Rabbit pAb
Remarks/Alias	EDRF; ERAF; AHSP
Species	Human
GeneID (human)	51327
GeneID	51327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-102 of human AHSP (NP_057717.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NZD4
Protein Weight	12kDa
Shipping	Ice bag