


Anti-Aif Polyclonal Antibody For Wb, If, Ihc, Elisa - 095831

Item Number: CM0016592



Introduction

This AIF polyclonal antibody (095831) is validated for WB, IF, IHC, ELISA. Reacts with human, mouse, rat AIF. Suitable for apoptosis and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	AIF Rabbit pAb
Remarks/Alias	AIF; AUNX1; CMT2D; CMTX4; COWCK; DFNX5; NADMR; NAMSD; PDCD8; COXPD6; SEMDHL
Species	Human
GeneID	9131
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 334-613 of human AIF (NP_004199.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095831
Protein Weight	67kDa
Shipping	Ice bag