


Anti-Aif Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O95831

Item Number: CM0003665



Introduction

AIF Rabbit Monoclonal Antibody validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat AIF (apoptosis-inducing factor 1, mitochondrial). Ideal for apoptosis and mitochondrial metabolism research.

[Learn More](#)

Parameter	Value
Product Name	AIF Rabbit mAb
Remarks/Alias	AIF; AUNX1; CMT2D; CMTX4; COWCK; DFNX5; NADMR; NAMSD; PDCD8; COXPD6; SEMDHL
Species	Human
Gene ID (Human)	9131
Gene ID	9131
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 450-550 of human AIF (O95831).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95831
Protein Weight	67kDa
Shipping	Ice bag