

Anti-Aif1/Iba1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P55008

Item Number: CM0007675



Introduction

Rabbit monoclonal antibody targeting human AIF1/IBA1 (SWISS P55008). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for macrophage activation, actin-binding, and inflammatory studies.

[Learn More](#)

Parameter	Value
Product Name	AIF1/IBA1 Rabbit PolymAb®
Remarks/Alias	IBA1; IRT1; AIF-1; IRT-1
Species	Human
GeneID (Human)	199
GeneID	199
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human AIF1/IBA1 (NP_001614.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55008
Protein Weight	17kDa
Shipping	Ice bag