

Anti-Aif1/Iba1 Mouse Monoclonal Antibody For Wb, If, Ihc-P, Elisa - O70200

Item Number: CM0012210



Introduction

Mouse monoclonal antibody against AIF1/IBA1 (AIF-1, Iba1) for WB, IF, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Detects ~17 kDa protein involved in actin binding and macrophage function. Ideal for inflammation and immunology research.

[Learn More](#)

Parameter	Value
Product Name	AIF1/IBA1 Mouse mAb
Remarks/Alias	G1; Iba1; AIF-1; D17H6S50E
Species	Mouse
Gene ID (Human)	199
Gene ID	11629
Immunogen	Recombinant protein of mouse AIF1/IBA1.
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O70200
Protein Weight	17kDa
Shipping	Ice bag